



## Intracellular Staining Using PFA and Saponin

### Reagents:

1X PBS supplemented with 2% FBS, 0.02%  $\text{NaN}_3$

10X PBS

Paraformaldehyde (PFA)

Saponin

Sodium carbonate

Distilled  $\text{H}_2\text{O}$

### Procedures:

1. Prepare PFA solution:
  - a. Heat up to boiling 90 mL of  $\text{H}_2\text{O}$  in a microwave.
  - b. Let temperature cool down.
  - c. Add 4g of PFA when temperature reaches  $68^\circ\text{C}$ .
  - d. Agitate with a magnetic stirrer at  $50^\circ\text{C}$ <sup>1</sup> and add some droplets of sodium carbonate<sup>2</sup>.
  - e. Allow the PFA to undergo complete dissolution (it becomes totally transparent).
  - f. Cool down the solution in ice.
  - g. Add 10 mL of 10x PBS.
  - h. Agitate with a magnetic stirrer until complete dissolution.
  - i. Adjust pH to 7,3.
  - j. Filter with a  $\varnothing 45 \mu\text{m}$  if necessary.
  - k. Aliquote in 5 mL volume and freeze at  $-20^\circ\text{C}$ <sup>3</sup>.
  - l. Dilute in 1x PBS to get a 2% concentration
  
2. Prepare PBS-saponin solution: 1X PBS, 5% FBS, 0,5% saponin. Dissolution of saponin should be done in agitation for 2 hours, in the dark, right before usage.

<sup>1</sup> Do not let temperature raise above  $65^\circ\text{C}$  to avoid PFA degradation.

<sup>2</sup> PFA can only undergo complete dissolution at basic pH.

<sup>3</sup> 4% PFA solution has a short shelf life. Use freshly prepared or thawed PFA solution for your experiments.

3. Transfer lymphocytes to a 15-mL falcon tube or to a 96-well plate.
4. Sediment lymphocytes by centrifugation at 1400 rpm for 5 min and discard the supernatant.
5. Perform surface staining as usual or proceed for the next step.

*Fixation step*

6. Resuspend the pellet in 1 mL (tubes) or 150  $\mu$ L (plate) 4% PFA.
7. Agitate for 5 minutes at room temperature.
8. Add 10 mL (tubes) or 150  $\mu$ L (plate) PBS and sediment cells by centrifugation at 1400 rpm for 5 min.
9. Discard the supernatant and repeat step 5.

*Permeabilization step*

10. Add 5 mL (tubes) or 300  $\mu$ L (plate) PBS-saponin.
11. Agitate for 10 minutes at room temperature.
12. Sediment cells by centrifugation at 1400 rpm for 5 min.
13. Add the intracellular antibodies diluted in PBS-saponin
14. Incubate for 30 minutes at 4°C, in the dark.
15. Add 5 mL (tubes) or 300  $\mu$ L (plate) PBS-saponin and sediment cells by centrifugation at 1400 rpm for 5 min.
16. Discard the supernatant and repeat step 5.
17. Discard the supernatant and resuspend cells in PBS.